

NAVIGATOR KIT

(COMMITTED TO EXCELLENCE - WORKING TOGETHER)

NATURAL LABOUR AND BIRTH KIT



**100%
EO
Sterilization**



NAVIGATOR KITS & DRAPRES HEALTHCARE PVT LTD

NO-2/2, VENKATRAMAN ROAD, CHINNA CHOKKIKULAAM, MADURAI 625002
PHONE NO:- 0452 2533669, 2533698, 4367889 MOBILE NO:- 98941-25965

NAVIGATOR NORMAL DELIVERY KIT - 100% EO STERILIZATION

The kit consists of all most all the vital items, essential for the Normal Labour & Spontaneous Delivery
All the product contents in the kit, are of superior quality procured from Multinational companies

KIT CONTENTS

MOTHER PART - 1 (PRE - DELIVERY)

SL	PRODUCT DESCRIPTION	QTY
1	DISPOSABLE MOTHER GOWN	1 NOS
2	ABSORBENT UNDER PAD	1 NOS
3	MATERNITY DRAPE SHEET	1 NOS
4	ABSORBENT UNDER PAD WITH FENESTRATION	1 NOS
5	CENTER HOLE SHEET 120 CMS X 80 CMS	1 NOS
6	LEGGINGS FOR MOTHER	1 PAIR
7	GAUZE SWABS 10 CMS X 10 CMS 12 PLY [10'S PACK]	1 PACK

A

SURGEON - PREPERATION PART

SL	PRODUCT DESCRIPTION	QTY
1	REINFORCED SURGEON GOWN	1 NOS
2	ABSORBENT TOWEL 40CMS X 40 CMS	1 NOS
3	SURGEON CAP	1 NOS
4	APRON FOR NURSE	1 NOS
5	NURSE CAP	1 NOS
6	FACE MASK - SURGEON (1 PCS) & NURSE (1 PCS)	2 PCS

B

BABY PART

SL	PRODUCT DESCRIPTION	QTY
1	TOWEL FOR BABY 80 CMS X 80 CMS	1 NOS
2	BABY BLANKET 1 MTR X 1 MTR	1 NOS
3	BABY WEAR	1 NOS

C

PROCEDURE PART

SL	PRODUCT DESCRIPTION	QTY
1	INSTRUCTIONAL SHEET	1 NOS
2	DIMENSIONAL SHEET	1 NOS
3	MEASURING TAPE	1 NOS
4	DISPOSABLE SAFETY RAZOR FOR MOTHER	1 NOS
5	STERILE GLOVES	4 PAIR
6	SCALPEL (STAINLESS STEEL) WITH PLASTIC HANDLE	1 NOS
7	UMBILICAL CORD CLAMP	1 NOS
8	INFANT MUCUS EXTRACTOR	1 NOS
9	DISPOSABLE SYRINGE WITH NEEDLE 10 ML	1 NOS
10	POVIDONE IODINE 100 ML BOTTLE	1 NOS
11	BIO-HAZARD - WASTE COLLECTION BAG	1 NOS

D

MOTHER PART - 2 (POST DELIVERY)

SL	PRODUCT DESCRIPTION	QTY
1	DISPOSABLE PAD FIXATOR [MESH PANTY]	1 NOS
2	ABSORBENT MEDI PAD [MATERNITY PAD]	2 NOS

E