



# Sterling

## Diagnositics & Therapy Equipments Pvt Ltd

F2-B Raksha Towers, No. 600/601 P.H. Road, Arumbakkam, Chennai - 106. India

Phone : +91 44 4218 7450, Fax : +91 44 4218 7313

E-mail : [info@sterlingdiagnostics.in](mailto:info@sterlingdiagnostics.in), Website : [www.sterlingdiagnostics.in](http://www.sterlingdiagnostics.in)

### **Patient Health Monitor**

#### **Features:**

- ✳ **12.1" inch color high definition (800 x 480) TFT screen**
- ✳ **Parameters: SPO2, ECG, NIBP, TEMP, RESP, ETCO2(optional)**
- ✳ **Wide-angle alarm indicator makes patients safer**
- ✳ **Easy operation of soft and colorful buttons**
- ✳ **Advanced patient information management system**
- ✳ **Connect with central monitoring system for remote and wireless monitor**
- ✳ **High capacity storage, power-off protection**
- ✳ **3-level visual/audio alarms and alarm events can be stored and reviewed**
- ✳ **Built-in rechargeable Lithium battery**
- ✳ **Equipped with internal printer**

