



# Sterling

## Diagnosics & Therapy Equipments Pvt Ltd

F2-B Raksha Towers, No. 600/601 P.H. Road, Arumbakkam, Chennai - 106. India

Phone : +91 44 4218 7450, Fax : +91 44 4218 7313

E-mail : [info@sterlingdiagnostics.in](mailto:info@sterlingdiagnostics.in), Website : [www.sterlingdiagnostics.in](http://www.sterlingdiagnostics.in)

### **Snore Stopper**

*Snore Stopper™ is a revolutionary device designed to reduce snoring through a patent nerve stimulation process. It features a specially designed biosensor with logic-control IC. When Snore Stopper™ detects snoring, a gentle, 4-second electrical pulse will be sent automatically to stimulate nerves on the snorer's skin. Then the brains senses this nerve stimulation, and remind or train the snorer him/herself to change the sleeping posture, thereby reducing snoring.*

*Snore Stopper™ features these advantages with absolutely no side effects!*

- ✳ *Smooth pulse without disturbing normal sleep*
- ✳ *Automatically detects snores*
- ✳ *Not affected if any buttons pressed while sleeping*
- ✳ *Automatic power-off after 12 hours of activation*
- ✳ *Fashionable watch-shaped design*

