



Sterling

Diagnosics & Therapy Equipments Pvt Ltd

F2-B Raksha Towers, No. 600/601 P.H. Road, Arumbakkam, Chennai - 106. India

Phone : +91 44 4218 7450, Fax : +91 44 4218 7313

E-mail : info@sterlingdiagnostics.in, Website : www.sterlingdiagnostics.in

Vital Signs Health Monitor

Features:

- ✳ **7 inch high resolution color TFT LCD screen**
- ✳ **Light weight(only 2.05Kg) and very portable**
- ✳ **Parameters: ECG,Resp,NIBP,Temp,SpO2,Pulse rate**
- ✳ **Adjustable 3-level audible and visual alarms**
- ✳ **Suitable for adult,pediatric and neonate**
- ✳ **Accurate SpO2 measurement in low perfusion situation**
- ✳ **Fast and accurate blood pressure in reading, reduce the problems of motion in NIBP measure**
- ✳ **Build-in WAN technology provides the ability to network with central monitoring for remote monitoring and diagnosis**

