

**PHILIPS OXYGEN CONCENTRATOR 10L**



## Specifications & Part Numbers

### Specifications

<b>Flow Rate:</b>	10 LPM in 1 liter increments
<b>Electrical Requirements:</b>	120V/60Hz ( $\pm 10\%$ )
<b>Oxygen Percentage Indicator:</b>	Green Light greater than 82% Yellow Light between 70-82% Red Light less than 70%
<b>Oxygen Concentration:</b>	92 $\pm$ 4% @ 8-10 LPM 94 $\pm$ 2% @ 3-7 LPM 92 $\pm$ 4% @ 1-2 LPM
<b>Weight:</b>	53 $\pm$ 0.5 lbs
<b>Dimensions:</b>	27" x 19" x 13"
<b>Storage/transport temperature:</b>	-30 to 160°F
<b>Operating temperature:</b>	55 to 90°F
<b>Storage/transport humidity:</b>	Up to 95%, non-condensing
<b>Operating pressure:</b>	10-30 psig
<b>Warranty:</b>	One year standard Optional two years (excludes compressor) Optional two year compressor
<b>Alarm indicators:</b>	High system pressure Low system pressure Power failure Low oxygen level (OPI@model available on the M10605 model only) No oxygen flow